



ADVALENE

Polypropylene

Monofilament Synthetic Non-Absorbable Surgical Suture

Reliable performance for confidence in contaminated wound closure





























ADVA TRAY™ - the confidence to suture with ease

ADVALENE Polypropylene

Monofilament Synthetic Non-Absorbable Surgical Suture

With exceptional strength comes confidence that lasts

Product Description

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA805	BLUE	0	70	1/2 CIRCLE		25	-	12
2	ADVAP824	BLUE	0	100	1/2 CIRCLE		36	-	12
3	ADVA830	BLUE	0	70	1/2 CIRCLE		30	HEAVY	12
4	ADVA833	BLUE	0	150	1/2 CIRCLE		40	-	12
5	ADVAP833	BLUE	0	90	1/2 CIRCLE		50	-	12
6	ADVA842	BLUE	0	70	1/2 CIRCLE		30	-	12
7	ADVA846	BLUE	0	90	1/2 CIRCLE		30	-	12
8	ADVA894	BLUE	0	90	1/2 CIRCLE		30	HEAVY	12
9	ADVA8014	BLUE	0	75	1/2 CIRCLE		26	-	12
10	ADVA8040	BLUE	0	70	1/2 CIRCLE		36	-	12
11	ADVA8041	BLUE	0	70	1/2 CIRCLE		40	-	12
12	ADVA8074	BLUE	0	100	1/2 CIRCLE		40	HEAVY	12
13	ADVA8075	BLUE	0	100	1/2 CIRCLE		40	-	12
14	ADVA831	BLUE	1	100	1/2 CIRCLE		60	-	12
15	ADVA834	BLUE	1	150	1/2 CIRCLE		48	HEAVY	12
16	ADVA834H	BLUE	1	150	1/2 CIRCLE		40	HEAVY	12
17	ADVAP834	BLUE	1	100	1/2 CIRCLE		50	-	12
18	ADVA840	BLUE	1	70	1/2 CIRCLE		30	HEAVY	12
19	ADVA843	BLUE	1	70	1/2 CIRCLE		40	HEAVY	12
20	ADVA883	BLUE	1	100	1/2 CIRCLE		45	HEAVY	12
21	ADVA896	BLUE	1	100	1/2 CIRCLE		30	HEAVY	12
22	ADVA8035	BLUE	1	100	1/2 CIRCLE		40	-	12
23	ADVA8038	BLUE	1	90	1/2 CIRCLE		75	-	12
24	ADVA8039	BLUE	1	70	1/2 CIRCLE		40	-	12
25	ADVA8044	BLUE	1	75	3/8 CIRCLE		40	-	12
26	ADVA8045	BLUE	1	90	1/2 CIRCLE		45	-	12
































 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

LP - Loop
































Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
27	ADVA8046	BLUE	1	100	1/2 CIRCLE		50	-	12
27	ADVA8118	BLUE	1	90	1/2 CIRCLE		30	-	12
29	ADVA8126	BLUE	1	70	3/8 CIRCLE		12	-	12
30	ADVA8129	BLUE	1	75	1/2 CIRCLE		40	-	12
31	ADVA8054	BLUE	1	75	1/2 CIRCLE		37	-	12
32	ADVA8073	BLUE	1	150	1/2 CIRCLE		40	-	12
33	ADVA8083	BLUE	1	150	1/2 CIRCLE		48	HEAVY	12
34	ADVA8084	BLUE	1	75	3/8 CIRCLE		36	-	12
35	ADVA8091	BLUE	1	90	1/2 CIRCLE		45	-	12
36	ADVA8116	BLUE	1	90	1/2 CIRCLE		24	-	12
37	ADVA8090	BLUE	2	76	1/2 CIRCLE		50	-	12
38	ADVA8130	BLUE	2	45	1/2 CIRCLE		26	-	12
39	ADVA806	BLUE	2-0	90	1/2 CIRCLE		26	-	12
40	ADVA807	BLUE	2-0	70	1/2 CIRCLE		26	-	12
41	ADVA819	BLUE	2-0	45	3/8 CIRCLE		15	-	12
42	ADVA841	BLUE	2-0	70	1/2 CIRCLE		30	-	12
43	ADVA844	BLUE	2-0	90	1/2 CIRCLE		30	-	12
44	ADVA850	BLUE	2-0	90	1/2 CIRCLE		36	-	12
45	ADVA863	BLUE	2-0	70	3/8 CIRCLE		45	-	12
46	ADVA8005	BLUE	2-0	90	1/2 CIRCLE		26	-	12
47	ADVA8008	BLUE	2-0	70	3/8 CIRCLE		30	-	12
48	ADVA8009	BLUE	2-0	75	1/2 CIRCLE		30	-	12
49	ADVA8026T	BLUE	2-0	45	3/8 CIRCLE		26	PRIME	24
50	ADVA8051	BLUE	2-0	45	3/8 CIRCLE		24	-	12
51	ADVA8076	BLUE	2-0	75	3/8 CIRCLE		20	-	12
52	ADVA8101	BLUE	2-0	75	3/8 CIRCLE		30	-	12
53	ADVA8111	BLUE	2-0	45	3/8 CIRCLE		26	PRIME	12
54	ADVA8115	BLUE	2-0	90	1/2 CIRCLE		24	-	12
55	ADVA8119	BLUE	2-0	90	1/2 CIRCLE		25	-	12
56	ADVA8121	BLUE	2-0	90	1/2 CIRCLE		30	-	12
57	ADVA8245	BLUE	2-0	45	1/2 CIRCLE		25	-	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut
































Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
58	ADVA8245L	BLUE	2-0	100	1/2 CIRCLE		25	-	12
59	ADVA8411	BLUE	2-0	75	1/2 CIRCLE		26	-	12
60	ADVA8631	BLUE	2-0	75	STRAIGHT		60	-	12
61	ADVA018	BLUE	3-0	70	3/8 CIRCLE		22	-	12
62	ADVA019	BLUE	3-0	70	3/8 CIRCLE		16	-	12
63	ADVAP422	BLUE	3-0	60	3/8 CIRCLE		26	-	12
64	ADVA800	BLUE	3-0	90	3/8 CIRCLE		25	PRIME	12
65	ADVA802	BLUE	3-0	70	3/8 CIRCLE		25	-	12
66	ADVA820	BLUE	3-0	45	3/8 CIRCLE		25	PRIME	12
67	ADVA821	BLUE	3-0	60	3/8 CIRCLE		25	-	12
68	ADVA822L	BLUE	3-0	75	1/2 CIRCLE		15	-	12
69	ADVA825	BLUE	3-0	90	1/2 CIRCLE		25	-	12
70	ADVA838	BLUE	3-0	70	1/2 CIRCLE		25	-	12
71	ADVA860	BLUE	3-0	75	STRAIGHT		60	-	12
72	ADVA887	BLUE	3-0	70	1/2 CIRCLE		20	-	12
73	ADVA8001	BLUE	3-0	70	1/2 CIRCLE		22	-	12
74	ADVA8003	BLUE	3-0	70	1/2 CIRCLE		36	-	12
75	ADVA8004	BLUE	3-0	70	1/2 CIRCLE		40	-	12
76	ADVA8006	BLUE	3-0	45	3/8 CIRCLE		26	PRIME	12
77	ADVA8021	BLUE	3-0	75	1/2 CIRCLE		22	-	12
78	ADVA8021T	BLUE	3-0	45	3/8 CIRCLE		26	P Curved	24
79	ADVA8042	BLUE	3-0	70	1/2 CIRCLE		24	-	12
80	ADVA8052	BLUE	3-0	75	3/8 CIRCLE		24	-	12
81	ADVA8061	BLUE	3-0	75	1/2 CIRCLE		16	-	12
82	ADVA8077	BLUE	3-0	75	3/8 CIRCLE		24	-	12
83	ADVA8078	BLUE	3-0	75	1/2 CIRCLE		18	-	12
84	ADVA8110	BLUE	3-0	45	3/8 CIRCLE		26	P Curved	12
85	ADVA8132	BLUE	3-0	70	3/8 CIRCLE		12	-	12
86	ADVA8832	BLUE	3-0	75	1/2 CIRCLE		26	LONG TAPER	12
87	ADVA815	BLUE	4-0	45	3/8 CIRCLE		16	-	12
88	ADVA816	BLUE	4-0	60	3/8 CIRCLE		16	PRIME	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut
































Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
89	ADVA848	BLUE	4-0	90	1/2 CIRCLE		16	LONG TAPER	12
90	ADVA849	BLUE	4-0	70	1/2 CIRCLE		16	CV	12
91	ADVA851	BLUE	4-0	70	3/8 CIRCLE		16	PRIME	12
92	ADVA857	BLUE	4-0	90	1/2 CIRCLE		17	-	12
93	ADVA870	BLUE	4-0	70	3/8 CIRCLE		16	PRIME	12
94	ADVA879	BLUE	4-0	45	3/8 CIRCLE		12	-	12
95	ADVA899	BLUE	4-0	70	1/2 CIRCLE		17	-	12
96	ADVA8007	BLUE	4-0	45	3/8 CIRCLE		19	PRIME	12
97	ADVA8043	BLUE	4-0	75	1/2 CIRCLE		26	-	12
98	ADVA8047	BLUE	4-0	60	3/8 CIRCLE		16	PRIME	12
99	ADVA8053	BLUE	4-0	75	3/8 CIRCLE		19	-	12
100	ADVA8058	BLUE	4-0	90	1/2 CIRCLE		13	-	12
101	ADVA8064	BLUE	4-0	45	3/8 CIRCLE		25	PRIME	12
102	ADVA8079	BLUE	4-0	150	3/8 CIRCLE		19	-	12
103	ADVA8092	BLUE	4-0	90	3/8 CIRCLE		26	-	12
104	ADVA8103	BLUE	4-0	60	3/8 CIRCLE		10	-	12
105	ADVA8112	BLUE	4-0	90	1/2 CIRCLE		19	-	12
106	ADVA8120	BLUE	4-0	75	1/2 CIRCLE		30	-	12
107	ADVA8122	BLUE	4-0	90	1/2 CIRCLE		30	-	12
108	ADVA8127	BLUE	4-0	70	3/8 CIRCLE		12	-	12
109	ADVA826	BLUE	5-0	70	3/8 CIRCLE		15	LONG TAPER	12
110	ADVA862	BLUE	5-0	70	3/8 CIRCLE		13	PRIME	12
111	ADVA880	BLUE	5-0	70	3/8 CIRCLE		12	PRIME	12
112	ADVA881	BLUE	5-0	70	3/8 CIRCLE		16	CV	12
113	ADVA881SP	BLUE	5-0	75	3/8 CIRCLE		18	-	12
114	ADVA889	BLUE	5-0	90	3/8 CIRCLE		12	PRIME	12
115	ADVA8010	BLUE	5-0	75	1/2 CIRCLE		18	-	12
116	ADVA8031	BLUE	5-0	75	1/2 CIRCLE		13	320 MICRONS	12
117	ADVA8048	BLUE	5-0	60	3/8 CIRCLE		16	PRIME	12
118	ADVA8055	BLUE	5-0	90	1/2 CIRCLE		17	PRIME	12
119	ADVA8065	BLUE	5-0	45	3/8 CIRCLE		20	PRIME	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut








Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
120	ADVA8080	BLUE	5-0	75	3/8 CIRCLE		12	-	12
121	ADVA8113	BLUE	5-0	90	1/2 CIRCLE		19	-	12
122	ADVA8123	BLUE	5-0	90	1/2 CIRCLE		30	-	12
123	ADVA8527	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	12
124	ADVA8872T	BLUE	5-0	45	3/8 CIRCLE		13	PRIME	24
125	ADVA8882	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	12
126	ADVA8882T	BLUE	5-0	45	3/8 CIRCLE		16	PRIME	24
127	ADVA861	BLUE	6-0	75	3/8 CIRCLE		8	-	12
128	ADVA868	BLUE	6-0	75	3/8 CIRCLE		12	-	12
129	ADVA878	BLUE	6-0	45	3/8 CIRCLE		12	-	12
130	ADVA878SP	BLUE	6-0	45	3/8 CIRCLE		15	-	12
131	ADVAP875	BLUE	6-0	75	1/2 CIRCLE		8	-	12
132	ADVA8060	BLUE	6-0	60	3/8 CIRCLE		10	-	12
133	ADVA8066	BLUE	6-0	75	3/8 CIRCLE		11	PRIME	12
134	ADVA8098	BLUE	6-0	60	3/8 CIRCLE		13	-	12
135	ADVA8114	BLUE	6-0	90	1/2 CIRCLE		19	-	12
136	ADVA8117	BLUE	6-0	45	3/8 CIRCLE		16	PRIME	12
137	ADVA8718	BLUE	6-0	60	3/8 CIRCLE		13	282 MICRONS	12
138	ADVA8868T	BLUE	6-0	45	3/8 CIRCLE		11	PRIME	24
139	ADVA8115A	BLUE	7-0	60	3/8 CIRCLE		9.5	-	12
140	ADVAP837	BLUE	7-0	75	1/2 CIRCLE		13	-	12
141	ADVA2076SA	BLUE	7-0	75	3/8 CIRCLE		6	-	12
142	ADVA2078L	BLUE	7-0	75	1/2 CIRCLE		8	-	12
143	ADVA8068	BLUE	7-0	60	3/8 CIRCLE		8	-	12
144	ADVA8068H	BLUE	7-0	60	3/8 CIRCLE		8	-	36
145	ADVA8069	BLUE	7-0	75	3/8 CIRCLE		8	-	12
146	ADVA8107	BLUE	7-0	90	1/2 CIRCLE		13	-	12
147	ADVA900SP	BLUE	8-0	38	3/8 CIRCLE		6	-	12
148	ADVA8033	BLUE	8-0	13	3/8 CIRCLE		5	-	12
149	ADVA8070	BLUE	8-0	60	3/8 CIRCLE		8	-	12
150	ADVA8102	BLUE	8-0	60	1/4 CIRCLE		6	-	12



















































 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
151	ADVA8131	BLUE	8-0	45	1/2 CIRCLE		8	-	12
152	ADVA2780	BLUE	9-0	13	3/8 CIRCLE		5	-	12
153	ADVA3709SP	BLUE	9-0	38	3/8 CIRCLE		6	MICRO POINT	12
154	ADVA8022	BLUE	9-0	23	3/8 CIRCLE		6.5	-	12
155	ADVA2790	BLUE	10-0	13	3/8 CIRCLE		4	-	12
156	ADVA2794	BLUE	10-0	13	3/8 CIRCLE		5	-	12
157	ADVA8049	BLUE	10-0	25	3/8 CIRCLE		6	-	12

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
158	ADVA824	BLUE	0	90	1/2 CIRCLE	 LT	 LT	26	12
159	ADVA8056	BLUE	0	75	1/2 CIRCLE			13	12
160	ADVA808	BLUE	2-0	70	1/2 CIRCLE			25	12
161	ADVA832	BLUE	2-0	70	STRAIGHT			30	12
162	ADVA837	BLUE	2-0	90	1/2 CIRCLE	 CV	 CV	17	12
163	ADVA877	BLUE	2-0	90	1/2 CIRCLE			25	12
164	ADVA898	BLUE	2-0	70	1/2 CIRCLE			17	12
165	ADVA8128	BLUE	2-0	70	1/2 CIRCLE			20	12
166	ADVA8526	BLUE	2-0	90	1/2 CIRCLE	 LT	 LT	31	12
167	ADVA8876	BLUE	2-0	90	1/2 CIRCLE			24	12
168	ADVA8886	BLUE	2-0	90	1/2 CIRCLE			26	12
169	ADVA8922	BLUE	2-0	90	3/8 CIRCLE			22	12
170	ADVA8937	BLUE	2-0	90	1/2 CIRCLE	 CV	 CV	17	12
171	ADVA8977	BLUE	2-0	90	1/2 CIRCLE			26	12
172	ADVA818	BLUE	3-0	60	3/8 CIRCLE	 LT	 LT	13	12
173	ADVA822	BLUE	3-0	70	1/2 CIRCLE	 LT	 LT	25	12
174	ADVA836	BLUE	3-0	90	1/2 CIRCLE			17	12
175	ADVA852	BLUE	3-0	90	1/2 CIRCLE	 LT	 LT	25	12
176	ADVA897	BLUE	3-0	70	1/2 CIRCLE			17	12
177	ADVA8002	BLUE	3-0	70	3/8 CIRCLE			17	12
178	ADVA8018	BLUE	3-0	90	3/8 CIRCLE			31	12
179	ADVA8036	BLUE	3-0	90	1/2 CIRCLE			25	12
180	ADVA8057	BLUE	3-0	75	1/2 CIRCLE			13	12
181	ADVA8085	BLUE	3-0	75	1/2 CIRCLE			18	12
182	ADVA8100	BLUE	3-0	90	1/2 CIRCLE			26	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
183	ADVA8108	BLUE	3-0	90	1/2 CIRCLE			25	12
184	ADVA8124	BLUE	3-0	70	1/2 CIRCLE			25	12
185	ADVA8332	BLUE	3-0	90	1/2 CIRCLE			30	12
186	ADVA8422	BLUE	3-0	90	1/2 CIRCLE			24	12
187	ADVA8522	BLUE	3-0	90	1/2 CIRCLE			26	12
188	ADVA8525	BLUE	3-0	90	1/2 CIRCLE	LT	LT	31	12
189	ADVA8558	BLUE	3-0	90	1/2 CIRCLE			17	12
190	ADVA8667	BLUE	3-0	90	1/2 CIRCLE			16	12
191	ADVA8936	BLUE	3-0	90	1/2 CIRCLE	CV	CV	17	12
192	ADVA8976	BLUE	3-0	90	1/2 CIRCLE			26	12
193	ADVA835	BLUE	4-0	70	1/2 CIRCLE	CV	CV	17	12
194	ADVA845	BLUE	4-0	90	1/2 CIRCLE	CV	CV	25	12
195	ADVA853	BLUE	4-0	90	1/2 CIRCLE			17	12
196	ADVA867	BLUE	4-0	90	1/2 CIRCLE			25	12
197	ADVA870DA	BLUE	4-0	70	3/8 CIRCLE	P	P	16	12
198	ADVA8015	BLUE	4-0	90	1/2 CIRCLE			22	12
199	ADVA8019	BLUE	4-0	90	1/2 CIRCLE			40	12
200	ADVA8020	BLUE	4-0	90	1/2 CIRCLE			36	12
201	ADVA8023	BLUE	4-0	90	1/2 CIRCLE			20	12
202	ADVA8026	BLUE	4-0	90	1/2 CIRCLE			16	12
203	ADVA8034	BLUE	4-0	90	1/2 CIRCLE			16	12
204	ADVA8086	BLUE	4-0	60	1/2 CIRCLE			13	12
205	ADVA8094	BLUE	4-0	90	1/2 & 3/8 CIRCLE	ADVABLACK	ADVABLACK	17	12
206	ADVA8099	BLUE	4-0	75	3/8 CIRCLE			8	12
207	ADVA8125	BLUE	4-0	70	1/2 CIRCLE			25	12
208	ADVA8134	BLUE	4-0	90	3/8 CIRCLE			20	12
209	ADVA8204	BLUE	4-0	60	1/2 CIRCLE	LT	LT	13	12
210	ADVA8204PL3	BLUE	4-0	60	1/2 CIRCLE	CV	CV	13	12
211	ADVA8354	BLUE	4-0	70	3/8 CIRCLE			16	12
212	ADVA8521	BLUE	4-0	90	1/2 CIRCLE	LT	LT	26	12
213	ADVA8521T	BLUE	4-0	90	1/2 CIRCLE	LT	LT	26	12
214	ADVA8557	BLUE	4-0	90	1/2 CIRCLE	LT	LT	17	12
215	ADVA8557PL3	BLUE	4-0	90	1/2 CIRCLE	CV LONG TAPER	CV LONG TAPER	17	12
216	ADVA8665	BLUE	4-0	75	1/2 CIRCLE			13	12

Round Bodied

Reverse Cutting

Cutting

Taper Cut











































































Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
217	ADVA8761	BLUE	4-0	90	1/2 CIRCLE	LT	LT	20	12
218	ADVA8845	BLUE	4-0	90	1/2 CIRCLE			26	12
219	ADVA8884T	BLUE	4-0	45	1/2 CIRCLE			16 & 26	24
220	ADVA8935	BLUE	4-0	90	1/2 CIRCLE			17	12
221	ADVA8935PL3	BLUE	4-0	90	1/2 CIRCLE	CV	CV	17	12
222	ADVA810	BLUE	5-0	75	1/2 CIRCLE	LT	LT	13	12
223	ADVA827	BLUE	5-0	70	3/8 CIRCLE	LT	LT	12	12
224	ADVA856	BLUE	5-0	90	1/2 CIRCLE	LT	LT	16	12
225	ADVA882	BLUE	5-0	70	3/8 CIRCLE	LT	LT	16	12
226	ADVA892	BLUE	5-0	70	3/8 CIRCLE	LT	LT	20	12
227	ADVA8016	BLUE	5-0	90	3/8 CIRCLE			17	12
228	ADVA8024	BLUE	5-0	90	1/2 CIRCLE			16	12
229	ADVA8027	BLUE	5-0	75	1/2 CIRCLE			13	12
230	ADVA8050	BLUE	5-0	90	1/2 CIRCLE			20	12
231	ADVA8096	BLUE	5-0	90	1/2 CIRCLE			18	12
232	ADVA8133	BLUE	5-0	70	3/8 CIRCLE			12	12
233	ADVA8205	BLUE	5-0	60	1/2 CIRCLE	LT	LT	13	12
234	ADVA8205PL3	BLUE	5-0	60	1/2 CIRCLE	CV	CV	13	12
235	ADVA8321	BLUE	5-0	60	3/8 CIRCLE	ADVABLACK	ADVABLACK	13	12
236	ADVA8556	BLUE	5-0	90	1/2 CIRCLE			17	12
237	ADVA8556PL3	BLUE	5-0	90	1/2 CIRCLE	CV	CV	17	12
238	ADVA8664	BLUE	5-0	60	1/2 CIRCLE			13	12
239	ADVA8710	BLUE	5-0	75	1/2 CIRCLE			13	12
240	ADVA8710PL3	BLUE	5-0	75	1/2 CIRCLE	CV	CV	13	12
241	ADVA8710S	BLUE	5-0	45	1/2 CIRCLE	LT	LT	13	12
242	ADVA8721	BLUE	5-0	75	3/8 CIRCLE			13	12
243	ADVA8803	BLUE	5-0	60	1/2 CIRCLE			11	12
244	ADVA8803PL3	BLUE	5-0	60	1/2 CIRCLE	CV	CV	11	12
245	ADVA8810	BLUE	5-0	60	1/2 CIRCLE			10	12
246	ADVA8011	BLUE	6-0	75	3/8 CIRCLE			13	12
247	ADVA8025	BLUE	6-0	60	1/2 CIRCLE			13	12
248	ADVA8028	BLUE	6-0	60	1/2 CIRCLE			8	12
249	ADVA8030	BLUE	6-0	60	1/2 CIRCLE			11	12
250	ADVA8032	BLUE	6-0	45	1/2 CIRCLE			8	12

Round Bodied

Reverse Cutting

Cutting

Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
251	ADVA8059	BLUE	6-0	60	3/8 CIRCLE			9.5	12
252	ADVA8062	BLUE	6-0	75	3/8 CIRCLE			10	12
253	ADVA8071	BLUE	6-0	75	3/8 CIRCLE			10	12
254	ADVA8081	BLUE	6-0	75	1/2 CIRCLE			10	12
255	ADVA8087	BLUE	6-0	75	3/8 CIRCLE			13	12
256	ADVA8104	BLUE	6-0	60	3/8 CIRCLE			13	12
257	ADVA8135	BLUE	6-0	90	3/8 CIRCLE			13	12
258	ADVA8807	BLUE	6-0	60	3/8 CIRCLE			9.5	12
259	ADVA8109	BLUE	6-0	60	3/8 CIRCLE			8	12
260	ADVA8206PL3	BLUE	6-0	60	1/2 CIRCLE			13	12
261	ADVA8207	BLUE	6-0	60	1/2 CIRCLE			10	12
262	ADVA8207PL3	BLUE	6-0	60	1/2 CIRCLE			10	12
263	ADVA8216	BLUE	6-0	75	1/2 CIRCLE			13	12
264	ADVA8307	BLUE	6-0	60	3/8 CIRCLE			13	12
265	ADVA8597	BLUE	6-0	60	3/8 CIRCLE			11	12
266	ADVA8681	BLUE	6-0	75	3/8 CIRCLE			13	12
267	ADVA8707	BLUE	6-0	60	3/8 CIRCLE			13	12
268	ADVA8717	BLUE	6-0	90	1/2 CIRCLE			17	12
269	ADVA8726	BLUE	6-0	60	3/8 CIRCLE			13	12
270	ADVA8726PL3	BLUE	6-0	60	3/8 CIRCLE			13	12
271	ADVA8802	BLUE	6-0	60	3/8 CIRCLE			11	12
272	ADVA8804PL3	BLUE	6-0	60	3/8 CIRCLE			10	12
273	ADVA811	BLUE	7-0	60	3/8 CIRCLE			6	12
274	ADVA2013	BLUE	7-0	75	3/8 CIRCLE			13	12
275	ADVA2070L	BLUE	7-0	75	3/8 CIRCLE			10	12
276	ADVA2078	BLUE	7-0	60	1/2 CIRCLE			10	12
277	ADVA2078B	BLUE	7-0	60	COMPOUND CURVE			8	12
278	ADVA2079L	BLUE	7-0	75	3/8 CIRCLE			9	12
279	ADVA8012	BLUE	7-0	60	3/8 CIRCLE			11	12
280	ADVA8017	BLUE	7-0	75	3/8 CIRCLE			9.5	12
281	ADVA8029	BLUE	7-0	60	1/2 CIRCLE			8	12
282	ADVA8037	BLUE	7-0	60	1/2 CIRCLE			8	12
283	ADVA8063	BLUE	7-0	60	3/8 CIRCLE			8	12
284	ADVA8067	BLUE	7-0	60	1/2 CIRCLE			9	12
285	ADVA8082	BLUE	7-0	75	1/2 CIRCLE			10	12
286	ADVA8088	BLUE	7-0	60	3/8 CIRCLE			9.5	12
287	ADVA8097	BLUE	7-0	75	3/8 CIRCLE			9	12

 Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut

Sr. No.	Double Arm Code	Suture			Needles (1 & 2)				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
288	ADVA8304	BLUE	7-0	60	3/8 CIRCLE	ADVABLACK	ADVABLACK	9.5	12
289	ADVA8702	BLUE	7-0	60	3/8 CIRCLE			9	12
290	ADVA8704	BLUE	7-0	60	3/8 CIRCLE			9.5	12
291	ADVA8725	BLUE	7-0	60	3/8 CIRCLE			13	12
292	ADVA8735	BLUE	7-0	60	3/8 CIRCLE			8	12
293	ADVA8735H	BLUE	7-0	60	3/8 CIRCLE			8	36
294	ADVA8747	BLUE	7-0	60	3/8 CIRCLE			9.5	12
295	ADVA8747H	BLUE	7-0	60	3/8 CIRCLE			9.5	36
296	ADVA2088L	BLUE	8-0	75	3/8 CIRCLE			8	12
297	ADVA8013	BLUE	8-0	60	3/8 CIRCLE			6	12
298	ADVA8072	BLUE	8-0	60	3/8 CIRCLE			9.5	12
299	ADVA8089	BLUE	8-0	60	3/8 CIRCLE			8	12
300	ADVA8093	BLUE	8-0	90	3/8 CIRCLE	ADVABLACK	ADVABLACK	8	12
301	ADVA8105	BLUE	8-0	45	3/8 CIRCLE	ADVAPASS	ADVAPASS	6.5	12
302	ADVA8106	BLUE	8-0	60	3/8 CIRCLE	ADVABLACK	ADVABLACK	8	12
303	ADVA8136	BLUE	8-0	90	1/2 CIRCLE			9	12
304	ADVA8632	BLUE	8-0	75	1/2 CIRCLE			6	12
305	ADVA8703	BLUE	8-0	60	3/8 CIRCLE	ADVAPASS	ADVAPASS	9	12
306	ADVA8732	BLUE	8-0	60	3/8 CIRCLE			6.5	12
307	ADVA8741	BLUE	8-0	60	3/8 CIRCLE	MICRO VASCULAR	MICRO VASCULAR	8	12
308	ADVA8741H	BLUE	8-0	60	3/8 CIRCLE	MICRO VASCULAR	MICRO VASCULAR	8	36
309	ADVA1710	BLUE	10-0	30	3/8 CIRCLE			6.5	12
310	ADVA1713	BLUE	10-0	23	STRAIGHT	MICROPOINT	MICROPOINT	16	12
311	ADVA8095	BLUE	10-0	75	STRAIGHT			16	12

⊙ Round Bodied

∇ Spatulated Micropoint

≡ Straight Cut

⊕ Taper Cut

△ Cutting

Product Description

Composition	Polypropylene
Coating	None
Sizes	U.S.P. Size 2 to 10-0
Colour	Blue
Absorption Profile	Permanent
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Taper Cut, Trocar, Reverse Cutting, Blunt Point, Oval Round Bodied, Diamond Point & Sternum Reverse Cutting

Monofilament Material



DURABLE COAT
RETAINS SHARPNESS...
PASS AFTER PASS

RIBBED BODY
FOR SUPERIOR CONTROL

CHAMFERING
MINIMUM SUTURE DETACHMENT


ELECTROPOLISHING
SMOOTH & BRIGHT

PREMIUM NEEDLES
ensure
UNSURPASSED SUTURING EXPERIENCE

Tensile Strength

86% more knot pull tensile strength Vs. U.S.P. limit ensures permanent wound support.

Knot pull tensile strength (KgF)¹



REFERENCE _____

¹ : Data on file with AMS - Bench test done on Polypropylene U.S.P. size 1, 40 mm Round bodied needles

The only Polypropylene with an advanced dispensing system

ADVATRAY™ - the confidence to
suture with ease

ADVATRAY™ a unique
dispensing system
designed to deliver ease
of suture dispensing
with minimal memory in
the suture strand⁷



Minimizing Suture Memory...
Maximizing Knot Security

**Enhanced robust packaging to ensure seal strength
integrity and sterility assurance throughout shelf life**

Advancing the Science of Cardiovascular Needles

Premium Needles for consistent Superior Performance

ADVAPASS TECHNOLOGY
ensures
UNSURPASSED SUTURING EXPERIENCE

ADVAPRIME SPECIALLY ENGINEERED DESIGN	Advanced point geometry and intelligent engineering assures steady, superior control and smoothness
ADVAGUARD BLUNT POINT	ADVAGUARD is engineered to diminish the chance of needle stick injuries, delivering both optimal strength and ductility with increased precision
ADVABLACK VISIBILITY ENABLING CONFIDENCE	ADVABLACK needles have been designed to offer outstanding visibility. Specially crafted for enhanced assurance and confidence in the operating zone

ADVAPASS advanced coating technology offers consistency in ease of penetration and maintains needle sharpness and performance from start to finish pass after pass

Minimal tissue trauma enables faster recovery



Advanced MedTech Solutions

P-21-22, 25-26 & 34-35, GIDC, Manjusr, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India
T : (91) 2667 671302 | amsltd.com | info@amsltd.com