



## ADVAMRYL

Poliglecaprone - 25 (Undyed & Violet)

Monofilament Synthetic Absorbable Surgical Suture

Supreme handling performance for confidence  
in subcuticular skin closure



**ADVATRAY™** - the confidence to  
suture with ease

ADVAMRYL

Poliglecaprone - 25 (Undyed & Violet)  
Monofilament Synthetic Absorbable Surgical Suture



Predictable strength for superior support

Product Description












































Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
1	ADVA1011	UNDYED	0	70	3/8 CIRCLE		30	-	12
2	ADVA1014	UNDYED	0	90	1/2 CIRCLE		48	-	12
3	ADVA1015	VIOLET	0	70	1/2 CIRCLE		37	-	12
4	ADVA1026	UNDYED	0	70	1/2 CIRCLE		40	HEAVY	12
5	ADVA1027	VIOLET	0	76	1/2 CIRCLE		35	-	12
6	ADVA1029	VIOLET	0	70	1/2 CIRCLE		37	-	12
7	ADVA1030	VIOLET	0	90	1/2 CIRCLE		40	HEAVY	12
8	ADVA1041	VIOLET	0	70	1/2 CIRCLE		40	HEAVY	12
9	ADVA1051	VIOLET	0	75	3/8 CIRCLE		20	-	12
10	ADVA1052	VIOLET	0	70	1/2 CIRCLE		30	-	12
11	ADVA1652	UNDYED	0	70	1/2 CIRCLE		25	-	12
12	ADVA1642	UNDYED	0	70	1/2 CIRCLE		30	-	12
13	ADVA1742	VIOLET	0	70	1/2 CIRCLE		30	-	12
14	ADVA3442	VIOLET	0	70	1/2 CIRCLE		31	LONG TAPER	12
15	ADVA3443	VIOLET	0	70	1/2 CIRCLE		36	LONG TAPER	12
16	ADVA3601	VIOLET	0	70	1/2 CIRCLE		36	-	12
17	ADVA1006	UNDYED	1	90	1/2 CIRCLE		40	-	12
18	ADVA1016	VIOLET	1	70	1/2 CIRCLE		37	-	12
19	ADVA1017	UNDYED	1	90	1/2 CIRCLE		48	-	12
20	ADVA1031	VIOLET	1	70	1/2 CIRCLE		40	HEAVY	12
21	ADVA1042	VIOLET	1	90	1/2 CIRCLE		40	HEAVY	12
22	ADVA1043	VIOLET	1	90	1/2 CIRCLE		50	BLUNT POINT	12
23	ADVA3457	VIOLET	1	70	1/2 CIRCLE		36	LONG TAPER	12
24	ADVA3755	VIOLET	1	70	1/2 CIRCLE		31	LONG TAPER	12
25	ADVA3759U	UNDYED	1	70	1/2 CIRCLE		40	-	12
26	ADVA1005	UNDYED	2-0	70	3/8 CIRCLE		30	-	12
27	ADVA1018	VIOLET	2-0	70	3/8 CIRCLE		24	-	12
28	ADVA1019	VIOLET	2-0	70	1/2 CIRCLE		26	LONG TAPER	12
29	ADVA1020	VIOLET	2-0	70	1/2 CIRCLE		37	-	12
30	ADVA1044	VIOLET	2-0	70	1/2 CIRCLE		35	-	12
31	ADVA1053	VIOLET	2-0	70	1/2 CIRCLE		30	-	12
32	ADVA1603	UNDYED	2-0	70	1/2 CIRCLE		25	-	12
33	ADVA1602	UNDYED	2-0	70	5/8 CIRCLE		26	-	12
34	ADVA1641	UNDYED	2-0	70	1/2 CIRCLE		30	-	12
35	ADVA1665	UNDYED	2-0	70	1/2 CIRCLE		26	-	12

Round Bodied

Oval Round Bodied

Reverse Cutting

Taper Cut

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
36	ADVA1666	UNDYED	2-0	70	1/2 CIRCLE		36	-	12
37	ADVA1666S	UNDYED	2-0	45	1/2 CIRCLE		36	-	12
38	ADVA1741	VIOLET	2-0	70	1/2 CIRCLE		30	-	12
39	ADVA3200	VIOLET	2-0	70	1/2 CIRCLE		31	LONG TAPER	12
40	ADVA3327	UNDYD	2-0	70	3/8 CIRCLE		26	-	12
41	ADVA3431	UNDYD	2-0	70	3/8 CIRCLE		26	-	12
42	ADVA3441	VIOLET	2-0	70	1/2 CIRCLE		31	-	12
43	ADVA945P	UNDYED	2-0	90	1/2 CIRCLE		36	-	12
44	ADVA1001	UNDYED	3-0	70	1/2 CIRCLE		20	-	12
45	ADVA1004	UNDYED	3-0	70	3/8 CIRCLE		24	-	12
46	ADVA1008	UNDYED	3-0	90	3/8 CIRCLE		19	-	12
47	ADVA1010	UNDYED	3-0	75	3/8 CIRCLE		26	-	12
48	ADVA1022	VIOLET	3-0	70	3/8 CIRCLE		24	-	12
49	ADVA1023	VIOLET	3-0	70	1/2 CIRCLE		30	-	12
50	ADVA1028	UNDYED	3-0	70	1/2 CIRCLE		23	-	12
51	ADVA1032	UNDYED	3-0	70	3/8 CIRCLE		24	-	12
52	ADVA1033	VIOLET	3-0	70	1/2 CIRCLE		26	LONG TAPER	12
53	ADVA1035	VIOLET	3-0	70	1/2 CIRCLE		22	LONG TAPER	12
54	ADVA1046	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
55	ADVA1049	VIOLET	3-0	90	3/8 CIRCLE		19	PRIME	12
56	ADVA1050	VIOLET	3-0	75	STRAIGHT		60	PRIME	12
57	ADVA1054	UNDYED	3-0	70	3/8 CIRCLE		20	PRIME	12
58	ADVA1021SA	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
59	ADVA1326	UNDYED	3-0	70	3/8 CIRCLE		25	PRIME	12
60	ADVA1326S	UNDYED	3-0	45	3/8 CIRCLE		25	PRIME	12
61	ADVA1326MR	VIOLET	3-0	70	3/8 CIRCLE		25	PRIME	12
62	ADVA1663	UNDYED	3-0	70	1/2 CIRCLE		22	-	12
63	ADVA1664	UNDYED	3-0	70	1/2 CIRCLE		26	-	12
64	ADVA1737	VIOLET	3-0	70	1/2 CIRCLE		20	-	12
65	ADVA3202	UNDYED	3-0	70	3/8 CIRCLE		19	-	12
66	ADVA3326M	UNDYED	3-0	70	3/8 CIRCLE		26	-	12
67	ADVA3326MS	VIOLET	3-0	45	3/8 CIRCLE		25	-	12
68	ADVA3665	VIOLET	3-0	75	1/2 CIRCLE		31	LONG TAPER	12
69	ADVA1003	UNDYED	4-0	70	3/8 CIRCLE		24	-	12
70	ADVA1007	UNDYED	4-0	75	1/2 CIRCLE		17	-	12
71	ADVA1009	UNDYED	4-0	45	3/8 CIRCLE		19	PRIME	12
72	ADVA1024	VIOLET	4-0	70	3/8 CIRCLE		19	-	12
73	ADVA1025	VIOLET	4-0	70	1/2 CIRCLE		17	-	12
74	ADVA1034	VIOLET	4-0	70	1/2 CIRCLE		22	LONG TAPER	12
75	ADVA1036	UNDYED	4-0	70	3/8 CIRCLE		19	-	12
76	ADVA1038	VIOLET	4-0	75	1/2 CIRCLE		20	LONG TAPER	12
77	ADVA1048	UNDYED	4-0	70	3/8 CIRCLE		19	PRIME	12
78	ADVA1205	UNDYED	4-0	70	3/8 CIRCLE		16	PRIME	12

 Round Bodied

 Oval Round Bodied

 Reverse Cutting

 Cutting

 Taper Cut



Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type	Length (mm)	Speciality	
79	ADVA1205MR	VIOLET	4-0	70	3/8 CIRCLE		16	PRIME	12
80	ADVA1648	UNDYED	4-0	70	3/8 CIRCLE		16	LONG TAPER	12
81	ADVA1748	VIOLET	4-0	70	3/8 CIRCLE		16	LONG TAPER	12
82	ADVA3201	UNDYED	4-0	45	3/8 CIRCLE		19	-	12
83	ADVA3205	UNDYED	4-0	45	3/8 CIRCLE		16	PRIME	12
84	ADVA845G	UNDYED	4-0	45	3/8 CIRCLE		16	PRIME	12
85	ADVA1002	UNDYED	5-0	45	3/8 CIRCLE		19	PRIME	12
86	ADVA1037	VIOLET	5-0	70	1/2 CIRCLE		17	-	12
87	ADVA1039	VIOLET	5-0	75	3/8 CIRCLE		13	LONG TAPER	12
88	ADVA1055	VIOLET	5-0	40	CURVED		20	-	12
89	ADVA3203	UNDYED	5-0	45	3/8 CIRCLE		13	PRIME	12
90	ADVA3204	UNDYED	5-0	45	3/8 CIRCLE		16	PRIME	12
91	ADVA3209	UNDYED	5-0	70	3/8 CIRCLE		19	-	12
92	ADVA493G	UNDYED	5-0	45	3/8 CIRCLE		13	PRIME	12
93	ADVA844G	UNDYED	5-0	45	3/8 CIRCLE		16	PRIME	12
94	ADVA1040	VIOLET	6-0	75	1/2 CIRCLE		13	LONG TAPER	12
95	ADVA3214	UNDYED	6-0	45	3/8 CIRCLE		11	PRIME	12
96	ADVA1047	VIOLET	2-0	150	-	-	-	-	12

Sr. No.	Code	Suture			Needle				Units Per Box
		Colour	Size	Length (cm)	Curvature	Type 1	Type 2	Length (mm)	
97	ADVA1012	VIOLET	0	75	1/2 CIRCLE			30	12
98	ADVA1021	VIOLET	3-0	70	1/2 CIRCLE & 3/8 CIRCLE			26	12


Round Bodied	Reverse Cutting	Cutting
--------------	-----------------	---------

Product Description	
Composition	Glycolide and ε-Caprolactone
Coating	None
Size	Size 1 to 6-0
Colour	Undyed and Violet
Absorption Profile	Essentially complete up to 90 days
Needle Availability	A wide variety of Round Bodied (Taper Point), Ball Point, Cutting, Trocar, Reverse Cutting, Taper Cut, Blunt Point, Oval Round Bodied & Diamond Point needles

Tensile Strength		
Days of Implantation	Tensile Strength Retention	
	DYED	UNDYED
7 Days (1 to 5-0)	60%	55%
7 Days (6-0)	60%	50%
14 Days (1 to 3-0 & 5-0)	40%	20%
14 Days (4-0 & 6-0)	40%	15%

66% more knot pull tensile strength Vs. recommended limit ensures excellent support throughout the critical wound healing period


#### Knot pull tensile strength (KgF)<sup>1</sup>




REFERENCE \_\_\_\_\_

1 : Data on file with AMS - Bench test done on Poliglecaprone - 25, size 3-0,

#### Monofilament Material



Violet



Undyed



Provides excellent support for a range of critical wound healing period to secure closure in a variety of tissues, procedures and patients.



P-21-22, 25-26 & 34-35, GIDC, Manjusar, Tal - Savli, Dist. Vadodara - 391775, Gujarat, India  
T : (91) 2667 671302 | [amsltd.com](http://amsltd.com) | [info@amsltd.com](mailto:info@amsltd.com)