


*Gliding the Science of
Flow Restoration*



ADVA GLIDE

PTCA BALLOON CATHETER

ADVANCED
INTERVENTIONAL
CARDIOLOGY
SOLUTIONS



BAR Ø	1.25 mm	1.50 mm	2.00 mm	2.25 mm	2.50 mm	2.75 mm	3.00 mm	3.25 mm	3.50 mm	4.00 mm	4.50 mm
4	1.15	1.48	1.92	2.20	2.46	2.68	2.92	3.17	3.47	3.92	4.38
6	1.24	1.52	2.05	2.28	2.54	2.76	3.04	3.27	3.54	4.04	4.50
8	1.26	1.58	2.11	2.31	2.59	2.82	3.09	3.34	3.66	4.14	4.61
10	1.27	1.64	2.14	2.37	2.66	2.88	3.15	3.40	3.71	4.22	4.69
12	1.28	1.74	2.18	2.41	2.74	2.94	3.21	3.45	3.78	4.29	4.75
14	1.30	1.78	2.22	2.45	2.82	2.99	3.25	3.49	3.88	4.36	4.80
16	1.31	1.85	2.26	2.49	2.86	3.03	3.29	3.53	3.97	4.44	4.86
18	1.33	2.00	2.30	2.53	2.90	3.07	3.34	3.57	4.04	4.54	4.93
20	1.34	2.10	2.38	2.58	2.98	3.12	3.40	3.62	4.15	4.65	-

Ordering Information Chart

Balloon Diameter (mm)	BALLOON LENGTH (mm)							
	10	12	15	20	25	30	35	40
1.25	AS-125-10	AS-125-12	AS-125-15	-	-	-	-	-
1.50	AS-150-10	AS-150-12	AS-150-15	AS-150-20	-	-	-	-
2.00	AS-200-10	AS-200-12	AS-200-15	AS-200-20	AS-200-25	-	-	-
2.25	AS-225-10	AS-225-12	AS-225-15	AS-225-20	AS-225-25	AS-225-30	-	-
2.50	AS-250-10	AS-250-12	AS-250-15	AS-250-20	AS-250-25	AS-250-30	AS-250-35	AS-250-40
2.75	AS-275-10	AS-275-12	AS-275-15	AS-275-20	AS-275-25	AS-275-30	AS-275-35	AS-275-40
3.00	AS-300-10	AS-300-12	AS-300-15	AS-300-20	AS-300-25	AS-300-30	AS-300-35	AS-300-40
3.25	AS-325-10	AS-325-12	AS-325-15	AS-325-20	AS-325-25	AS-325-30	AS-325-35	AS-325-40
3.50	AS-350-10	AS-350-12	AS-350-15	AS-350-20	AS-350-25	AS-350-30	AS-350-35	AS-350-40
4.00	AS-400-10	AS-400-12	AS-400-15	AS-400-20	AS-400-25	AS-400-30	AS-400-35	AS-400-40
4.50	-	AS-450-12	-	AS-450-20	AS-450-25	AS-450-30	AS-450-35	AS-450-40

