



## ADVAKIT C-Section

Supreme handling performance for complete ease and confidence in Caesarean surgeries



Kit Code	Brand Name	Code	Suture		Needles	Foil Per Kit	Kit Per Box
			Size	Length (cm)			
AG01	ADVACRYL	ADVA2347DA180	1	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVAMRYL	ADVA1326	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME	01	
AG02	ADVACRYL	ADVA2347DA180	1	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVACRYL	ADVA2317	2-0	90	30 mm, 1/2 Circle, Round Bodied	01	
	ADVACAT	ADVA4259	1	76	40 mm, 1/2 Circle, Taper Point, HEAVY	01	
AG03	ADVACRYL	ADVA2347DA180	1	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVAMRYL	ADVA1326S	3-0	45	25 mm, 3/8 Circle, Cutting, PRIME	01	
	ADVAPD	ADVA9352	1	90	Fascia Closure	01	
AG04	ADVACRYL	ADVA2347DA180	1	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVACRYL	ADVA2317	2-0	90	30 mm, 1/2 Circle, Round Bodied	01	
	ADVAMRYL	ADVA3326M	3-0	70	26 mm, 3/8 Circle, Reverse Cutting	01	
AG06	ADVACRYL	ADVA2347DA180	1	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVAMRYL	ADVA1326S	3-0	45	25 mm, 3/8 Circle, Cutting, PRIME	01	
AG07	ADVACRYL	ADVA2347	1	90	40 mm, 1/2 Circle, Round Bodied, HEAVY	01	06
	ADVACRYL	ADVA2346	0	90	40 mm, 1/2 Circle, Round Bodied	01	
	ADVAMRYL	ADVA1326	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME	01	
AG08	ADVACRYL	ADVA2347DA180	1	180	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVAMRYL	ADVA1326	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME	01	
	ADVACAT	ADVA4259	1	76	40 mm, 1/2 Circle, Taper Point, HEAVY	01	
AG09	ADVACRYL	ADVA2347DA200	1	200	40 mm, 1/2 Circle, Taper Point & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVAMRYL	ADVA1326S	3-0	45	40 mm, 1/2 Circle, Round Bodied	01	
AG10	ADVACRYL	ADVA2346MM	0	110	40 mm, 1/2 Circle, Round Bodied	01	06
	ADVACRYL	ADVA2317	2-0	90	30 mm, 1/2 Circle, Round Bodied	01	
	ADVACAT	ADVA4241	2-0	76	30 mm, 1/2 Circle, Taper Point	01	
AG11	ADVACRYL	ADVA2347DA140	1	140	40 mm, 1/2 Circle, Round Bodied & 40 mm, 1/2 Circle, Taper Cut, HEAVY	01	06
	ADVAMRYL	ADVA1326S	3-0	45	25 mm, 3/8 Circle, Cutting, PRIME	01	
AG12	ADVACRYL	ADVA2360L2DA	1	200	40 mm, 1/2 Circle, Taper Point	01	06
	ADVAMRYL	ADVA1326S	3-0	45	25 mm, 3/8 Circle, Cutting, PRIME	01	
HS01	ADVACRYL	ADVA2347	1	90	40 mm, 1/2 Circle, Round Bodied, HEAVY	05	06
	ADVACRYL	ADVA2338	0	90	30 mm, 1/2 Circle, Round Bodied	01	
	ADVAMRYL	ADVA1326	3-0	70	25 mm, 3/8 Circle, Cutting, PRIME	01	

